

## **Appendix 2 *ancestryinfer* pipeline**

Output of *mixnmatch* simulations *or* your own data can be input into the *ancestryinfer* pipeline to run local ancestry inference following (Corbett-Detig & Nielsen, 2017).

### **Install**

*Option 1 – install dependencies:*

```
git clone https://github.com/Schumerlab/ancestryinfer.git
```

To install dependencies, follow instructions outlined in:

```
installation_instructions.txt
```

*Option 2 – load docker file for dependencies:*

With docker:

```
git clone https://github.com/Schumerlab/ancestryinfer.git
```

```
docker pull schumer/mixnmatch-ancestryinfer-image:mixnmatch-ancestryinfer-docker
```

```
docker run -it mixnmatch-ancestryinfer-docker bash
```

With singularity:

```
git clone https://github.com/Schumerlab/ancestryinfer.git
```

```
singularity pull docker://schumer/mixnmatch-ancestryinfer-image:mixnmatch-ancestryinfer-docker
```

```
singularity run mixnmatch-ancestryinfer-image_mixnmatch-ancestryinfer-docker.sif bash
```

### **Setting parameters in the configuration file**

There are example configuration files available on github:

```
hmm_configuration_file_parallel.cfg
```

```
hmm_configuration_file_nonparallel.cfg
```

Parameter descriptions:

Parameter	Description	Example	Include if
-----------	-------------	---------	------------

<code>genome1=</code>	User provided fasta file for species 1	<code>genome1=xiphophorus_birchmanni_10x_12Sep2018_yDAA6.fasta</code>	Always
<code>genome2=</code>	User provided fasta file for species 2	<code>genome2=Xmalinche_dovetail_assembly.fa</code>	Always
<code>read_type=</code>	Indicated whether data is paired end or single end	<code>read_type=PE</code>	Always
<code>read_list=</code>	Provide list (including full paths) to the reads to be analyzed	<code>read_list=combined_all_call_hybrids_read_list</code>  example list format for paired end data (single end file should contain one line per individual):  <code>./reads/CALL1_read1.fq.gz ./reads/CALL1_read2.fq.gz</code> <code>./reads/CALL2_read1.fq.gz ./reads/CALL2_read2.fq.gz</code> <code>./reads/CALL3_read1.fq.gz ./reads/CALL3_read2.fq.gz</code>	Always
<code>read_length=</code>	Provide expected read length	<code>read_length=150</code>	Always
<code>prop_genome_genome1_parent=</code>	Expected proportion of the genome derived from the parent species listed under genome1	<code>prop_genome_genome1_parent=0.5</code>	If not provided, AncestryHMM will attempt to estimate (may increase run time)
<code>number_indiv_per_job=</code>	Parallelize jobs such that each job processes this number of individuals.  Low numbers mean high parallelization and high number mean low parallelization.	<code>number_indiv_per_job=1</code>	Always
<code>program_path=</code>	Path to the program install folder	<code>program_path=/home/groups/schumer/shared_bin/Ancestry_HMM_pipeline</code>	If not provided the program will assume necessary scripts and

			programs are in the working directory
<code>provide_AIMs=</code>	Coordinates and identities of ancestry informative sites that distinguish the two parent species	<code>provide_AIMs=Xbirmanni10xgenome_ancestry_informative_sites_filterF1</code>  Example list format:  <pre> ScyDAA6-2-HRSCAF-26  58345  T      C ScyDAA6-2-HRSCAF-26  58976  T      A ScyDAA6-2-HRSCAF-26  59896  T      C ScyDAA6-2-HRSCAF-26  60164  G      A ScyDAA6-2-HRSCAF-26  63105  G      A ScyDAA6-2-HRSCAF-26  65532  G      A ScyDAA6-2-HRSCAF-26  66290  C      A ScyDAA6-2-HRSCAF-26  68233  T      C ScyDAA6-2-HRSCAF-26  70398  G      A ScyDAA6-2-HRSCAF-26  73869  G      A </pre>	Required unless provided genomes are on the same coordinate system and can be auto detected
<code>provide_counts=</code>	Counts of parental allele frequencies at ancestry informative sites (and recombination rates between adjacent sites if available)	<code>provide_counts=Xbirmanni10xgenome_Xmalinche_observed_parental_counts_filterF1</code>  Example format:  <pre> ScyDAA6-2-HRSCAF-26  163722  129  3  0  54  0.00000078 ScyDAA6-2-HRSCAF-26  166158  135  5  0  54  0.00001374 ScyDAA6-2-HRSCAF-26  166535  6    0  0  6  0.00000754 </pre> Columns are:  <pre> Chromosome site allele1_count_parent1 allele2_count_parent1 allele1_count_parent2 allele2_count_parent2 recombination_rate </pre>	If not provided, the program will assume that provided ancestry informative sites are fixed between species
<code>per_site_error=</code>	Per-site error parameter for HMM (i.e. due to sequencing error, contamination, etc)	<code>per_site_error=0.02</code>	Always
<code>gen_initial_admix=</code>	Estimated generation of initial admixture	<code>gen_initial_admix=20</code>	If not provided, AncestryHMM will attempt to estimate (may increase run time)
<code>focal_chrom_list=</code>	Provide a list of chromosomes to run (other	<code>focal_chrom_list=mychrs.txt</code>  Example:	Not required

	chromosomes will not be run)	ScyDAA6-2-HRSCAF-26 ScyDAA6-7-HRSCAF-50	
rec_M_per_bp=	Estimated recombination rate in Morgans/bp	rec_M_per_bp=0.00000002	Always; Use an estimate for a related species if not available
max_alignments=	Limit analysis to a maximum number of alignments (for computational speed)	max_alignments=2000000	Optional
retain_intermediate_files=	Keep all intermediate files. Warning: setting this to 1 results in a high space footprint for a large run; only recommended for troubleshooting.	retain_intermediate_files=0	Options are 1 to keep or 0 to delete.
posterior_thresh=	Posterior probability threshold to use for identifying ancestry transition intervals	posterior_thresh=0.9	Recommended 0.8-1
job_submit_command=	Option to run sequentially if using Docker image for dependencies or from a desktop computer. Set bash to run sequentially and sbatch to run in parallel on a slurm cluster	job_submit_command=bash or job_submit_command=sbatch	Always required

<code>slurm_command_map=</code>  <code>slurm_command_variant_call=</code>  <code>slurm_command_hmm=</code>	If running on a slurm cluster, provide cluster specific parameters for queues, time & memory	<code>slurm_command_map=#!/bin/sh #SBATCH --ntasks=1 #SBATCH --cpus-per-task=1 #SBATCH -p schumer --mem=64000 #SBATCH --time=02:30:00</code>  <code>slurm_command_variant_call=#!/bin/sh #SBATCH --ntasks=1 #SBATCH --cpus-per-task=1 #SBATCH -p schumer --mem=64000 #SBATCH --time=05:00:00</code>  <code>slurm_command_hmm=#!/bin/sh #SBATCH --ntasks=1 #SBATCH --cpus-per-task=1 #SBATCH -p schumer --mem=64000 #SBATCH --time=03:00:00</code>	Required if running on a cluster
--	--	---	----------------------------------

## Examples

Several example files are available with the git repository including example configuration files

## Running the pipeline

After setting the parameters in the configuration file and loading required dependencies, simply run:

```
perl mixnmatch/ simulate_admixed_genomes_v6.pl
hybrid_simulation_configuration.cfg
```

where path is the path to your simulator install

## **Appendix 3: Protocol used to generate low-coverage sequenced data for F<sub>1</sub> and F<sub>2</sub> individuals (*X. birchmanni* x *X. malinche*)**

### **1. Pre-charge the Tn5 with the adaptors.**

\_\_\_\_\_ Combine on ice, *in the following order*, and mix by pipetting up and down after each addition:

15 µl Tn5 (100 ng / µl)  
122 µl reassociation buffer/Glycerol (1:1 mix of reassociation buffer and glycerol, made in advance)  
3 µl each adaptor (1 and 2)

\_\_\_\_\_ Incubate in a thermal cycler at 37°C for 30 minutes (with a heated lid).

*Note: this will exceed the allowed volume per well for most thermocyclers, so the mixture should be split between three wells with the appropriate volume per well*

\_\_\_\_\_ While the Tn5 is pre-charging, add 3 µl of DNA (3-10 ng /µl) to each well on a 96 well plate (label “tagmentation”)

### **2. Tagmentation**

\_\_\_\_\_ Make up MasterMix for tagmentation reaction on ice. The volumes per 96 samples (107µL/well mmix in a strip) are:

120 ul Precharged Tn5 from step 1  
240 ul 5X TAPS Buffer (contains DMF) – collect waste separately  
500 ul H2O

\_\_\_\_\_ Add 7 uL of the MasterMix to each well of a 96 well plate, already containing your DNA. Mix up and down as you add the MasterMix. Spin briefly.

\_\_\_\_\_ Incubate at 55°C for 7 minutes.

### **3. Kill the Tn5**

\_\_\_\_\_ Add 2.5 µl 0.2% SDS to each reaction, mix up and down, spin briefly.

\_\_\_\_\_ Incubate at 55°C for 7 minutes in a thermal cycler.

### **4. PCR**

\_\_\_\_\_ Transfer 3 µl of your digested DNA from each well of your old 96 well plate to a new 96 well plate