Appendix 2 ancestryinfer pipeline

Output of *mixnmatch* simulations <u>or</u> your own data can be input into the *ancestryinfer* pipeline to run local ancestry inference following (Corbett-Detig & Nielsen, 2017).

Install

cd ancestryinfer

```
Option 1 – install dependencies:
git clone https://github.com/Schumerlab/ancestryinfer.git
To install dependencies, follow instructions outlined in:
installation instructions.txt
Test that the install and pipeline are working:
cd ancestryinfer
perl Ancestry HMM parallel v5.pl
hmm configuration file nonparallel.cfg
Option 2 – load docker file for dependencies:
With docker:
docker pull schumer/mixnmatch-ancestryinfer-image:mixnmatch-
ancestryinfer-docker
docker run -it mixnmatch-ancestryinfer-image bash
With singularity:
singularity pull docker://schumer/mixnmatch-ancestryinfer-
image:mixnmatch-ancestryinfer-docker
git clone https://github.com/Schumerlab/ancestryinfer.git
singularity run mixnmatch-ancestryinfer-image mixnmatch-ancestryinfer-
docker.sif bash
Test that the install and pipeline are working:
```

```
perl Ancestry_HMM_parallel_v5.pl
hmm_configuration_file_nonparallel.cfg
```

Setting parameters in the configuration file

There is are example configuration files available on github:

 ${\tt hmm_configuration_file_parallel.cfg}$

hmm_configuration_file_nonparallel.cfg

Parameter descriptions:

Parameter	Description	Example	Include if
genome1=	User provided fasta file for species 1	<pre>genome1=xiphophorus_birchmanni_10x_12Sep2018_yD AA6.fasta</pre>	Always
genome2=	User provided fasta file for species 2	<pre>genome2=Xmalinche_dovetail_assembly.fa</pre>	Always
read_type=	Indicated whether data is paired end or single end	read_type=PE	Always
read_list=	Provide list (including full paths) to the reads to be analyzed	read_list=combined_all_call_hybrids_read_list example list format for paired end data (single end file should contain one line per individual): ./reads/CALL1_read1.fq.gz ./reads/CALL1_read2.fq.gz ./reads/CALL2_read1.fq.gz ./reads/CALL2_read2.fq.gz ./reads/CALL3_read1.fq.gz ./reads/CALL3_read2.fq.gz	Always
read_length=	Provide expected read length	read_length=150	Always
<pre>prop_genome_g enome1_parent =</pre>	Expected proportion of the genome derived from the parent species listed under genome1	<pre>prop_genome_genome1_parent=0.5</pre>	If not provided, AncestryH MM will attempt to estimate (may increase run time)
number_indiv_ per_job=	Parallelize jobs such that each	number_indiv_per_job=1	Always

	1	T	1
	job processes		
	this number of		
	individuals.		
	Low numbers		
	mean high		
	parallelization		
	and high number		
	mean low		
	parallelization.		
program_path=	Path to the	<pre>program_path=/home/groups/schumer/shared_bin/An</pre>	If not
1 0 1		cestry_HMM_pipeline	provided the
	program install		-
	folder		program
			will assume
			necessary
			scripts and
			programs
			are in the
			working
			_
provide_AIMs=	0 1 1	<pre>provide_AIMs=Xbirchmanni10xgenome_ancestry_info</pre>	directory
bi ovine_ATM2=	Coordinates and	rmative_sites_filterF1	Required
	identities of		unless
	ancestry	Example list format:	provided
	informative sites		genomes are
	that distinguish	ScyDAA6-2-HRSCAF-26 58345 T C	on the same
	the two parent	ScyDAA6-2-HRSCAF-26 58976 T A	coordinate
	species	ScyDAA6-2-HRSCAF-26 59896 T C	system and
	species	ScyDAA6-2-HRSCAF-26 60164 G A ScyDAA6-2-HRSCAF-26 63105 G A	can be auto
		ScyDAA6-2-HRSCAF-26 65532 G A	detected
		ScyDAA6-2-HRSCAF-26 66290 C A	detected
		ScyDAA6-2-HRSCAF-26 68233 T C	
		ScyDAA6-2-HRSCAF-26 70398 G A	
mmand da assert		ScyDAA6-2-HRSCAF-26 73869 G A	TC /
<pre>provide_count s=</pre>	Counts of	<pre>provide_counts=Xbirchmanni10xgenome_Xmalinche_o bserved parental counts filterF1</pre>	If not
	parental allele	03c. 1ca_par circu1_counc3_f11ccf11	provided,
	frequencies at	Example format:	the program
	ancestry	r	will assume
	informative sites	ScyDAA6-2-HRSCAF-26 163722 129 3 0 54 0.00000078	that
	(and	ScyDAA6-2-HRSCAF-26 166158 135 5 0 54 0.00001374	provided
	recombination	ScyDAA6-2-HRSCAF-26 166535 6 0 0 6 0.00000754	ancestry
	rates between	Columns are:	informative
		Columnia are.	
	adjacent sites if	Chromosome site allele1_count_parent1	sites are
	available)	allele2_count_parent1 allele1_count_parent2	fixed
		allele2_count_parent2 recombination_rate	between
			species
per_site_erro	Per-site error	per_site_error=0.02	Always
r=	parameter for		
	HMM (i.e. due		
	THYINT (I.C. duc		

	to sequencing error, contamination, etc)		
gen_initial_a dmix=	Estimated generation of initial admixture	<pre>gen_initial_admix=20</pre>	If not provided, AncestryH MM will attempt to estimate (may increase run time)
<pre>focal_chrom_l ist=</pre>	Provide a list of chromosomes to run (other chromosomes	focal_chrom_list=mychrs.txt Example:	Not required
	will not be run)	ScyDAA6-2-HRSCAF-26	
rec M per bp=	Estimated	ScyDAA6-7-HRSCAF-50 rec M per bp=0.00000002	Always;
	recombination	'	Use an
	rate in		estimate for
	Morgans/bp		a related
			species if
			not
may alignment	T : 1	may alignments 2000000	available
<pre>max_alignment s=</pre>	Limit analysis to a maximum	max_alignments=2000000	Optional
	number of		
	alignments (for		
	computational		
	speed)		
retain_interm	Keep all	retain_intermediate_files=0	Options are
ediate_files=	intermediate		1 to keep or
	files. Warning:		0 to delete.
	setting this to 1		
	results in a high space footprint		
	for a large run;		
	only		
	recommended		
	for		
	troubleshooting.		
<pre>posterior_thr esh=</pre>	Posterior	posterior_thresh=0.9	Recommend
C311-	probability		ed 0.8-1
	threshold to use		
	for identifying		

	ancestry transition intervals		
job_submit_co mmand=	Option to run sequentially if using Docker image for dependencies or from a desktop computer. Set bash to run sequentially and sbatch to run in parallel on a slurm cluter	<pre>job_submit_command=bash or job_submit_command=sbatch</pre>	Always required
slurm_command _map= slurm_command _variant_call = slurm_command _hmm=	If running on a slurm cluster, provide cluster specific parameters for queues, time & memory	<pre>slurm_command_map=#!/bin/sh #SBATCHntasks=1 #SBATCHcpus-per-task=1 #SBATCH -p schumer mem=64000 #SBATCHtime=02:30:00 slurm_command_variant_call=#!/bin/sh #SBATCH ntasks=1 #SBATCHcpus-per-task=1 #SBATCH -p schumermem=64000 #SBATCHtime=05:00:00 slurm_command_hmm=#!/bin/sh #SBATCHntasks=1 #SBATCHcpus-per-task=1 #SBATCH -p schumer mem=64000 #SBATCHtime=03:00:00</pre>	Required if running on a cluster

Examples

Several example files are available with the git repository including example configuration files

Running the pipeline

After setting the parameters in the configuration file and loading required dependencies, simply run:

```
perl mixnmatch/simulate_admixed_genomes_v6.pl
hybrid_simulation_configuration.cfg
```

where path is the path to your simulator install

Appendix 3: Protocol used to generate low-coverage sequenced data for F_1 and F_2 individuals (X. birchmanni x X. malinche)

1. Pre-charge the Tn5 with the adaptors.

Combine on ice, <i>in the following order</i> , and mix by pipetting up and down after each addition:
$15~\mu l$ Tn5 (100 ng / $\mu l)$ $122~\mu l$ reassociation buffer/Glycerol (1:1 mix of reassociation buffer and glycerol, made in advance) $3~\mu l$ each adaptor (1 and 2)
Incubate in a thermal cycler at 37°C for 30 minutes (with a heated lid).
Note: this will exceed the allowed volume per well for most thermocyclers, so the mixture should be split between three wells with the appropriate volume per well
While the Tn5 is pre-charging, add 3 μ l of DNA (3-10 ng $/\mu$ l) to each well on a 96 well plate (label "tagmentation")
2. <u>Tagmentation</u>
Make up MasterMix for tagmentation reaction on ice. The volumes per 96 samples (107uL/well mmix in a strip) are: 120 ul Precharged Tn5 from step 1 240 ul 5X TAPS Buffer (contains DMF) – collect waste separately 500 ul H20
Add 7 uL of the MasterMix to each well of a 96 well plate, already containing your DNA. Mix up and down as you add the MasterMix. Spin briefly.
Incubate at 55°C for 7 minutes.
3. Kill the Tn5
Add 2.5 μl 0.2% SDS to each reaction, mix up and down, spin briefly.
Incubate at 55°C for 7 minutes in a thermal cycler.
4. <u>PCR</u>
Transfer 3 μl of your digested DNA from each well of your old 96 well plate to a new 96 well plate